



# SAFESEPT<sup>MAX</sup> Surface Disinfection FV

Surface disinfection • virucidal

- alcohol based
- bactericidal incl. MRSA, yeasticidal
- virucidal incl. Polio-, Adeno-, Rota- and Noro-viruses



## PRODUCT DESCRIPTION:

SAFESEPT MAX® Surface Disinfection FV is a ready-to-use, alcohol based (ethanol) rapid disinfectant. The product is specifically designed for applications requiring a targeted virucidal disinfection of alcohol resistant surfaces. Thanks to its innovative SLT booster technology, SAFESEPT MAX® Surface Disinfection FV also ensures an efficient disinfection in cases of high organic soiling (dirty conditions).

## AREAS OF APPLICATION:

Acc. to MDD: Rapid disinfection of alcohol-resistant, non-invasive medical devices.

Acc. to Biocidal Products Regulations: Rapid cleaning and disinfection of alcohol-resistant medical inventory and all types of surfaces.

## SPECTRUM OF ACTIVITY AND DOSAGE:

Bactericidal incl. MRSA, tuberculocidal, yeasticidal and virucidal incl. Polio-, Adeno-, Rota- and Noro viruses. Surface disinfection (wipe disinfection) (bactericidal, yeasticidal, dirty conditions): 1 min.; Tuberculocidal (*M. terrae*) acc. to EN 14348: 30 sec.; Virucidal acc. EN 14476 incl. HBV/HIV/HCV: 1 min.; Polio- and Noro virus (MNV): 1 min.; Adeno- and Rota virus: 30 sec.

## APPLICATION:

SAFESEPT MAX® Surface Disinfection FV can be applied by wiping (the recommended method of application) or by spraying. Surface disinfection (wipe disinfection): Apply a sufficient quantity of SAFESEPT MAX® Surface Disinfection FV on a suitable single use cloth and wipe surfaces. Ensure that the surfaces remain completely wetted with the disinfectant for the required contact time. For spray-based disinfection (recommended only for hard-to-reach surfaces), spray SAFESEPT MAX® Surface Disinfection FV onto the surface from a distance of at least 30 cm. Ensure that the surfaces remain completely wetted for the required contact time. Remove heavy soiling prior to disinfection using a single use cloth soaked in disinfectant. Afterwards, apply SAFESEPT MAX® Surface Disinfection FV to the surfaces to be disinfected. Once the contact time has expired, the surfaces can be wiped dry using a clean disposable cloth.

## COMPOSITION:

100 g contains: 60 g ethanol (CAS: 64-17-5); contains nonionic surfactants < 5%, fragrances.



# SAFESEPT<sup>MAX</sup> Surface Disinfection FV


Surface disinfection • full virucidal

- alcohol based
- bactericidal incl. MRSA, yeasticidal
- virucidal incl. Polio-, Adeno-, Rota- and Noro-viruses



### SPECIAL INFORMATION:

Highly flammable liquid and vapour. Causes serious eye irritation. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. Do not breathe vapour. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. Dispose of contents/container to approved disposal company or local collection.



### AVAILABLE ACCESSORIES:

- REF 979-6277 Sprayhead for 1 Liter Bottle
- REF 979-6276 Tab/dosing device for 5L canister
- REF 979-6278 Canister-key
- REF 979-6272 Dosing pump for 5L Canister

### PRODUCT STATUS/LISTINGS:

- Conform to the EU Biocidal Product Regulation (BPR, 528/2012).
- Conform to the EU guideline 93/42/EC for medical devices.

### AVAILABLE DELIVERY FORMS:

- REF 988-2742 1L Spray bottle
- REF 988-2743 5L Canister