



EuroSept®Plus Thermo Rinse

- Rinse agent for mechanical reprocessing of medical instruments
- Optimal wetting of items to be washed
- Ensures spot-free drying
- Miele and Melag approved



PRODUCT DESCRIPTION:

Henry Schein selected a cleaner liquid, a neutraliser and a rinsing liquid to accompany the Alkine powder for the perfect mechanical instrument cleaning cycle for by using thermodisinfectors. The liquids need to be added to their respective containers for automatic dosing during the cycle. The powder has to be added manually. EuroSept®Plus products are accepted by the major manufacturers of the thermodisinfectors.

Concentrated rinsing agent for spot free dried instruments in the last cycle of the mechanical reprocessing of instruments.

CONCENTRATION:

0.02-0.06%; dosing is performed with the aid of the dosing device/window provided in the thermodisinfectior. The fill level must be checked regularly (see machine manufacturer's instructions for use)

COMPOSITION:

5-15% non-ionic surfactants, corrosion inhibitors, solubilisers

SPECIAL INFORMATION:

Can causes skin irritation and serious eye damage. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. Protect from frost, heat and sunlight. Store in original container at room temperature.

AVAILABLE DELIVERY FORMS:

REF 988-2721 1L bottle

ALSO AVAILABLE:

- REF 988-2722 EuroSept®Plus Thermo Powder
- REF 988-2723 EuroSept®Plus Thermo Neutralizer
- REF 988-2724 EuroSept®Plus Thermo Cleanser